

## **PRODUCT INFORMATION**

Common Name MEDI-6570; MEDI 6570; MEDI 6570

Conjugate Unconjugated

**Synonyms** LOX1;LOXIN;SLOX1;CLEC8A;SCARE1;OLR1

Applications N/A
Recommended N/A
Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

Host Species N/A
IgG type IgG1
Reactivity Human
Target LOX-1
Uniprot ID P78380

**Description** Anti-LOX-1(golocdacimab biosimilar ☐ lgG1) mAb

**Delivery** 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store

**Storage & Shipping** Intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

**Usage** Research use only

