

**PRODUCT INFORMATION**

|   |   |
|---|---|
| <b>Common Name</b>                      | JNJ-63733657;JNJ 63733657;JNJ63733657;JNJ-3657;JNJ3657;JNJ 3657   |
| <b>Conjugate</b>                        | Unconjugated  |
| <b>Synonyms</b>                         | TAU;FTD1;MSTD;PPND;DDPAC;MAPTL;MTBT1;MTBT2;tau-40;FTDP-17;PPP1R103;Tau-PHF6   |
| <b>Applications</b>                     | N/A   |
| <b>Recommended Dilutions</b>            | N/A   |
| <b>Formulation &amp; Reconstitution</b> | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.  |
| <b>Host Species</b>                     | N/A   |
| <b>IgG type</b>                         | IgG1  |
| <b>Reactivity</b>                       | Human   |
| <b>Target</b>                           | MAPT  |
| <b>Uniprot ID</b>                       | P10636  |
| <b>Description</b>                      | Anti-MAPT(posdinemab biosimilar) [IgG1] mAb   |
| <b>Delivery</b>                         | 4-6 weeks   |
| <b>Storage&amp;Shipping</b>             | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| <b>Background</b>                       | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.   |
| <b>Usage</b>                            | Research use only   |

DIMABIO CONFIDENTIAL

