

PRODUCT INFORMATION

Common Name	IPH-2201, IPH2201, NN-8765, NNC-0141-0000-0100, NNC-141-01000, anti-NKG2A, humZ270
Conjugate	Unconjugated
Synonyms	CD159A;NKG2;NKG2A
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Humanized
IgG type	Human IgG4 - Kappa
Reactivity	Human
Target	NKG2A
Uniprot ID	P26715
Description	Anti-NKG2A(monalizumab biosimilar) mAb
Delivery	In Stock
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.






Figure 1. Flow cytometry analysis with 1 μ g/mL Anti-NKG2A (monalizumab biosimilar) mAb (BME100085) on Expi293 cells transfected with Human NKG2A protein and Human CD94 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

