

PRODUCT INFORMATION

Common Name AMG-145; AMG145; Repatha

Conjugate Unconjugated

FH3;PC9;FHCL3;NARC1;LDLCQ1;NARC-1;HCHOLA3 **Synonyms**

Applications N/A Recommended N/A

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A IgG type IgG2 Reactivity Human **Target** PCSK9 **Uniprot ID** Q8NBP7

Description Anti-PCSK9(evolocumab biosimilar IgG2) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

