

## **PRODUCT INFORMATION**

**Common Name** 

ANB 011;ANB011;ANB-011;GSK-4057190;GSK4057190;JEMPERLI;TSR 042;TSR042;TSR-042;WBP 285;WBP285;WBP-285

Conjugate

**Synonyms** PD1;PD-1;CD279;SLEB2;hPD-1;hPD-l;hSLE1

**Applications** Recommended N/A

**Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of Formulation & Reconstitution

reconstitution.

**Host Species** IgG type lgG4 Reactivity Human **Target** PD-1 **Uniprot ID** Q15116

Description Anti-PD-1(dostarlimab biosimilar ☐ IgG4) mAb

**Delivery** 4-6 weeks

Storage & Shipping

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Research grade biosimilar. Not for use in therapeutic or

> Email: info@dimabio.com Website: www.dimabio.com

**Background** diagnostic procedures for humans or animals.

Usage Research use only

