

## **PRODUCT INFORMATION**

APL 501; APL-501; APL501; CBT **Common Name** 

501;CBT501;CBT-501;GB 226;GB-226;GB226

Conjugate Unconjugated

**Synonyms** PD1;PD-1;CD279;SLEB2;hPD-1;hPD-l;hSLE1

**Applications** N/A Recommended N/A

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

**Host Species** N/A IgG type IgG4 Reactivity Human PD-1 **Target Uniprot ID** Q15116

**Description** Anti-PD-1(geptanolimab biosimilar IgG4) mAb

**Delivery** 

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage** 

