

PRODUCT INFORMATION

Common Name LVGN-3616;LVGN3616;LVGN 3616

Conjugate Unconjugated

Synonyms PD1;PD-1;CD279;SLEB2;hPD-1;hPD-l;hSLE1

Applications N/A
Recommended N/A
Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

Host Species N/A
IgG type IgG4
Reactivity Human
Target PD-1
Uniprot ID Q15116

Description Anti-PD-1(pradusinstobart biosimilar ☐ lgG4) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store

Storage & Shipping Intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

