

PRODUCT INFORMATION

AEX-1188;AEX1188;AEX

1188;INCMGA-0012;INCMGA

Common Name 0012;INCMGA0012;INCMGA12;INCMGA-12;INCMGA

12;MGA-012;MGA012;MGA 012;MGA-12;MGA12;MGA 12;Zynyz

Conjugate Unconjugated

PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1 **Synonyms**

Applications Recommended N/A

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species N/A IgG type IgG4 Reactivity Human PD-1 **Target Uniprot ID** Q15116

Description Anti-PD-1(retifanlimab biosimilar IgG4) mAb

Delivery 4-6 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

