Cat. No. BME100461



PRODUCT INFORMATION

Common Name CK 301;CK301;CK-301;TG 1501;TG1501;TG-1501

Conjugate Unconjugated

Synonyms B7-H;B7H1;PDL1;PD-L1;hPD-L1;PDCD1L1;PDCD1LG1

Applications Recommended N/A

Dilutions

Formulation &

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for Reconstitution

specific instructions of reconstitution.

Host Species N/A IgG type lgG1 Reactivity Human PD-L1 **Target** Q9NZQ7 **Uniprot ID**

Description Anti-PD-L1(cosibelimab biosimilar IgG1) mAb

Delivery 4-6 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use

> > Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping within a month, aliquot and store at -80°C (Avoid

repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans or

animals.

Usage Research use only

