Cat. No. BME100679



## **PRODUCT INFORMATION**

**Common Name** HS636;HS-636;HS 636

Conjugate Unconjugated

**Synonyms** B7-H;B7H1;PDL1;PD-L1;hPD-L1;PDCD1L1;PDCD1LG1

**Applications** Recommended

N/A Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for Formulation & Reconstitution

specific instructions of reconstitution.

**Host Species** N/A IgG type lgG1 Reactivity Human PD-L1 **Target** Q9NZQ7 **Uniprot ID** 

Description Anti-PD-L1(sudubrilimab biosimilar IgG1) mAb

**Delivery** 4-6 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use

> > Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping within a month, aliquot and store at -80°C (Avoid

repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans or

animals.

Usage Research use only

