

PRODUCT INFORMATION

Cejemly;CS-1001;CS1001;CS

N/A

1001;Sugelimab;WBP-3155;WBP3155;WBP **Common Name**

3155;Zegema;Zegemet;Zegemex;ZEJIEMEI

Conjugate Unconjugated

Synonyms B7-H;B7H1;PDL1;PD-L1;hPD-L1;PDCD1L1;PDCD1LG1

Applications N/A

Recommended

Background

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for Formulation & Reconstitution specific instructions of reconstitution.

Host Species N/A lgG4 IgG type

Reactivity Human **Target** PD-L1 **Uniprot ID** Q9NZQ7

Description Anti-PD-L1(sugemalimab biosimilar IgG4) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized

form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid Storage & Shipping

repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or

animals.

Usage Research use only



