

PRODUCT INFORMATION

Clone ID	37A11
Target	Rabbit IgG
Synonyms	N/A
Host Species	Rabbit
Description	Anti-Rabbit IgG antibody(37A11), Goat mAb
Delivery	In Stock
Uniprot ID	N/A
IgG type	Goat IgG
Clonality	Monoclonal
Reactivity	Rabbit
Applications	ELISA
Recommended Dilutions	ELISA 1/5000-10000
Purification	Purified from cell culture supernatant by affinity chromatography
Endotoxin	Less than 1.0 EU/μg by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μm) prior to use.
Background	N/A
Usage	Research use only
Conjugate	Unconjugated
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scr

SKU	Clone ID	Species Reactivity			
		Rabbit IgG	Human IgG	Rat IgG	Mouse IgG
DME101138	34B1	+++	++	-	-
DME101139	37A11	++	-	-	-



Figure 1. ELISA examination of Goat anti-rabbit IgG mAbs binding to immunoglobulins of different species.

DIMABIO CONFIDENTIAL

