

PRODUCT INFORMATION

Common Name r-anti-D Conjugate Unconjugated

Applications Recommended Dilutions

Formulation & Reconstitution $Lyophilized from sterile \ PBS, pH\ 7.4.\ Normally\ 5\%-8\%\ trehalose\ is\ added\ as\ protectants\ before\ lyophilization.$ Please see Certificate of Analysis for specific instructions of reconstitution.

Host Species lgG1 IgG type Reactivity Human Target RHD Uniprot ID Q02161

Description Anti-RHD(trinbelimab biosimilar IgG1) mAb

Delivery 4-6 weeks

use within a are shipped at . . or animals: Storage & Shipping

Background

Usage

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

Email: info@dimabio.com Website: www.dimabio.com

