

PRODUCT INFORMATION

Common Name	UC-961,cirmtuzumab
Synonyms	ROR1;NTRKR1
Applications	ELISA; Flow Cyt
Recommended Dilutions	ELISA 1:5000-10000; Flow Cyt 1:100
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	ROR1
Uniprot ID	Q01973
Description	Anti-ROR1 (zilovetamab biosimilar) mAb
Delivery	In Stock
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only



Anti-ROR1 (zilovetamab biosimilar) mAb ELISA
0.1 µg of Human ROR1, His tagged protein per well




Figure 1. ELISA plate pre-coated by 2 µg/ml (100 µl/well) Human ROR1, His tagged protein ([getskuurl sku="PME100399"]) can bind Anti-ROR1 Neutralizing antibody (**BME100073**) in a linear range of 0.64-16 µg/ml.




Figure 2. Flow cytometry analysis with Anti-ROR1 (zilovetamab biosimilar) mAb 15 µg/mL on Expi293 cells transfected with Human ROR1 (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

