Cat. No. BME100681



PRODUCT INFORMATION

BGB-A425;BGB A425 **Common Name**

Conjugate Unconjugated

TIM3;CD366;KIM-3;SPTCL;TIMD3;Tim-3;TIMD-3;HAVcr-2 **Synonyms**

Applications N/A Recommended

Dilutions

N/A

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species lgG1 IgG type Reactivity Human TIM3 **Target Uniprot ID** Q8TDQ0

Description Anti-TIM3(surzebiclimab biosimilar IgG1) mAb

Delivery

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid

Storage & Shipping repeated freezing and thawing). Lyophilized proteins

are shipped at ambient temperature.

Research grade biosimilar. Not for use in therapeutic **Background**

or diagnostic procedures for humans or animals.

Email: info@dimabio.com Website: www.dimabio.com

Usage Research use only

