Cat. No. BME100724



PRODUCT INFORMATION

INCAGN2390;INCAGN02390;INCAGN-2390;INCAGN-02390 **Common Name**

Conjugate Unconjugated

TIM3:CD366:KIM-3:SPTCL:TIMD3:Tim-3:TIMD-3:HAVcr-2 **Synonyms**

Applications N/A Recommended N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions Formulation & Reconstitution

of reconstitution.

Host Species N/A IgG type lgG1 Reactivity Human TIM3 **Target Uniprot ID** Q8TDQ0

Description Anti-TIM3(verzistobart biosimilar [IgG1) mAb

Delivery 4-6 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated

> > Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping freezing and thawing). Lyophilized proteins are shipped

at ambient temperature.

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background**

Research use only Usage

