

PRODUCT INFORMATION

Common Name	GEN1029;GEN-1029;Hx-DR5-01;Hx-DR5-05;iDD 004;iDD-004;iDD004
Conjugate	Unconjugated
Synonyms	DR5;CD262;KILLER;TRICK2;TRICKB;ZTNFR9;TRAILR2;TRICK2A;TRICK2B;TRAIL-R2;KILLER/DR5
Applications	N/A
Endotoxin	Less than 1.0 EU/μg by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	IgG1
Reactivity	Human
Target	TNFRSF10B
Uniprot ID	O14763
Description	Anti-TNFRSF10B(benufutamab biosimilar) IgG1) mAb
Delivery	4-6 weeks
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μm) prior to use.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

DIMABIO CONFIDENTIAL

