

PRODUCT INFORMATION

Common Name TBC

Conjugate Unconjugated

p75;TBPII;TNFBR;TNFR2;CD120b;TNFR1B;TNFR80;TNF-R75;p75TNFR;TNF-R-II **Synonyms**

Applications Recommended Dilutions N/A

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species N/A IgG type lgG1 Reactivity Human Target TNFRSF1B **Uniprot ID** P20333

Anti-TNFRSF1B(emunkitug biosimilar∏lgG1) mAb Description

Delivery

Storage & Shipping

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Email: info@dimabio.com Website: www.dimabio.com

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background**

Research use only Usage

