Cat. No. BME100762



## **PRODUCT INFORMATION**

**Common Name** BCD-180;BCD180;BCD 180

Conjugate Unconjugated

TCRBV1S1A1N1;TCRBV9S1 **Synonyms** 

**Applications** N/A Recommended N/A **Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

**Host Species** N/A IgG type lgG1 Reactivity Human **Target** TRBV9

**Uniprot ID** A0A0B4J1U6

**Description** Anti-TRBV9(seniprutug biosimilar ☐ IgG1) mAb

**Delivery** 4-6 weeks

Storage & Shipping

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

Email: info@dimabio.com Website: www.dimabio.com

Lyophilized proteins are shipped at ambient


témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only



/+86-400-006-0995(China)

