

PRODUCT INFORMATION

Common Name VGL101;VGL-101 Conjugate Unconjugated

PLOSL2;TREM-2;Trem2a;Trem2b;Trem2c **Synonyms**

Applications N/A Recommended N/A **Dilutions**

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

Host Species N/A IgG type lgG1 Reactivity Human **Target** TREM2 **Uniprot ID** Q9NZC2

Description Anti-TREM2(iluzanebart biosimilar IgG1) mAb

Delivery

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store Storage & Shipping

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

