

PRODUCT INFORMATION

Common Name TBC

Conjugate Unconjugated

Synonyms EGP1;GP50;M1S1;EGP-1;TROP2;GA7331;GA733-1;TACSTD2

Applications N/A Recommended Dilutions

N/A

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions

of reconstitution.

Host Species N/A IgG type lgG1 Reactivity Human **Target** Trop2 P09758 **Uniprot ID**

Description Anti-Trop2(tizetatug biosimilar IgG1) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping

temperature.

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background**

Research use only Usage

