

PRODUCT INFORMATION

BSI-045B;BSI-04502;BSI045B;BSI04502;TQC2731;TQC-2731 Common Name

Conjugate Unconjugated

TSLP Synonyms Applications N/A Recommended N/A Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of Formulation & Reconstitution

reconstitution.

Host Species IgG type lgG1 Reactivity Human TSLP **Target Uniprot ID** 0969D9

Description Anti-TSLP(bosakitug biosimilar IgG1) mAb

Delivery 4-6 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and

> > Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping

thawing). Lyophilized proteins are shipped at ambient temperature.

Research grade biosimilar. Not for use in therapeutic or **Background**

diagnostic procedures for humans or animals.

Usage Research use only

