

## **PRODUCT INFORMATION**

APX003;APX-003;BD0801;BD-0801;TK001;TK-001;Sevacizumab **Common Name** 

(not WHO assigned INN)

Conjugate Unconjugated

VPF:VEGF:MVCD1:L-VEGF **Synonyms** 

N/A **Applications** Recommended Dilutions

N/A

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please Formulation &

Reconstitution see Certificate of Analysis for specific instructions of

reconstitution.

**Host Species** N/A IgG type lgG1 Reactivity Human **VEGFA Target Uniprot ID** P15692

**Description** Anti-VEGFA(suvemcitug biosimilar [IgG1) mAb

**Delivery** 4-6 weeks

Storage & Shipping

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Email: info@dimabio.com Website: www.dimabio.com

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background** 

Usage Research use only

