

PRODUCT INFORMATION

SSS-23;SSS23;SSS 23;TTAC-0001;TTAC0001;TTAC **Common Name**

0001;TTAC-1;TTAC1;TTAC 1

Conjugate Unconjugated

Synonyms FLK1;CD309;VEGFR;VEGFR2

Applications N/A Recommended N/A

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species N/A IgG type IgG1 Reactivity Human VEGFR2 **Target Uniprot ID** P35968

Description Anti-VEGFR2(olinvacimab biosimilar [IgG1) mAb

4-6 weeks **Delivery**

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**

