

PRODUCT INFORMATION

| | |
|-----------------------------|--|
| Target | CD46, CD73, CDH17, CEACAM5, EPCAM, FZD10, GPA33, GUCY2C, HER3 |
| Description | The embedded and sliced tumor sections (4µm section) from Balb/C nu mice implanted with HT55 cells |
| Host Cells | HT55 |
| Applications | IHC |
| Host mice | Balb/C nu |
| Delivery | 1 week |
| Guarantee | 36 months from the date of manufacture |
| Background | These carefully prepared tissue slides offer researchers standardized and high-quality samples for antibody testing. Utilizing tumor models derived from immunocompromised mice, researchers can confidently assess antibody specificity and sensitivity, driving advancements in cancer diagnostics and therapeutics. |
| Storage&Shipping | Room temperature |
| Usage | Research use only |

DIMABIO CONFIDENTIAL

