

PRODUCT INFORMATION

Common Name	MGA017;MGA-017;MGA17;MGA-17;MGA 017;MGA 17;MGA018;MGA-018;MGA18;MGA-18;MGA 018;MGA 18;AEX4089;AEX-4089;AEX 4089
Conjugate	Biotinylated
Synonyms	B7H3;B7-H3;B7RP-2;4Ig-B7-H3;CD276
Applications	N/A
Endotoxin	Less than 1.0 EU/μg by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	IgG1
Reactivity	Human
Target	B7-H3
Uniprot ID	Q5ZPR3
Description	Biotinylated Anti-B7-H3(vobramitamab biosimilar) mAb
Delivery	6-8 weeks
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

