

PRODUCT INFORMATION

| | |
|------------------------------|---|
| Common Name | ARGX-110;ARGX110;JNJ 74494550;JNJ74494550;JNJ-74494550;JNJ-4550;JNJ4550 |
| Conjugate | Biotinylated |
| Synonyms | CD27L;LPFS3;CD27-L;CD27LG;TNFSF7;TNLG8A |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | N/A |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | CD70 |
| Uniprot ID | P32970 |
| Description | Biotinylated Anti-CD70(cusatuzumab biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage & Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

