

PRODUCT INFORMATION

Common Name	DR-01;DR01;DR 01
Conjugate	Biotinylated
Synonyms	KP43;KLRD1
Applications	N/A
Recommended Dilutions	15 µg/ml
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	IgG1
Reactivity	Human
Target	CD94
Uniprot ID	Q13241
Description	Biotinylated Anti-CD94(dibotatug biosimilar) mAb
Delivery	In Stock
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.






Figure 1. Flow cytometry analysis with 15µg/ml Biotinylated Anti-CD94(dibotatug biosimilar) mAb on HEK293 cells transfected with Human CD94 protein (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

