

PRODUCT INFORMATION

|                              |   |
|------------------------------|---|
| Common Name                  | 5A11  |
| Conjugate                    | Biotinylated  |
| Synonyms                     | CD355   |
| Applications                 | Flow Cyt  |
| Recommended Dilutions        | Flow Cyt 1:100  |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.  |
| Host Species                 | Humanized   |
| IgG type                     | IgG1 (K97R,L117A,L118A,E239D,M241L)   |
| Reactivity                   | Human   |
| Target                       | CRTAM   |
| Uniprot ID                   | O95727  |
| Description                  | Biotinylated Anti-CRTAM(Oxford Biotherapeutics 5A11) mAb  |
| Delivery                     | 2-3 weeks   |
| Storage & Shipping           | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature. |
| Background                   | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.   |
| Usage                        | Research use only   |

