

PRODUCT INFORMATION

| | |
|------------------------------|---|
| Common Name | DS-1024;DS1024;S-95026;S95206;Sym-004 (Combination with Futuximab);Sym004 (Combination with Futuximab) |
| Conjugate | Biotinylated |
| Synonyms | ERBB;ERRP;HER1;mENA;ERBB1;NNCIS;PIG61;NISBD2 |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | N/A |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | EGFR |
| Uniprot ID | P00533 |
| Description | Biotinylated Anti-EGFR(modotuximab biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage & Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

