

PRODUCT INFORMATION

Common Name	DE-766;DE766;DE 766;KI-0501;KI0501;KI0501;KI-501;KI 501;KI501;KI-0502;KI0502;KI0502;KI-502;KI502;KI502;KI502;OSAG101;OSAG-101;OSAG101;Taixinsheng;TheraCIM;Theraloc;VECTHIX
Conjugate	Biotinylated
Synonyms	ERBB;ERRP;HER1;mENA;ERBB1;NNCIS;PIG61;NISBD2
Applications	N/A
Endotoxin	Less than 1.0 EU/μg by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	IgG1
Reactivity	Human
Target	EGFR
Uniprot ID	P00533
Description	Biotinylated Anti-EGFR(nimotuzumab biosimilar) mAb
Delivery	6-8 weeks
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μm) prior to use.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

