Delivery

Uniprot ID

Background



PRODUCT INFORMATION

Clone ID 1A6

Target Exatecan

Synonyms N/A
Host Species Rabbit

DescriptionBiotinylated Anti-Exatecan antibody(1A6), Rabbit

mAb In Stock N/A

IgG type Rabbit IgG
Clonality Monoclonal
Reactivity Human
Applications ELISA

Recommended Dilutions ELISA 1/5000-10000

PurificationPurified from cell culture supernatant by affinity

chromatography

Formulation & Powder Reconstitution

Storage & Shipping

Storage & Shipping

Storage & Shipping

Storage & Shipping

intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

Exatecan (Exa) is a very potent inhibitor of topoisomerase I and anticancer agent. It has been intensively studied as a single agent, a large macromolecular conjugate and as the

Email: info@dimabio.com Website: www.dimabio.com

large macromolecular conjugate and as the payload component of antigen-dependent

antibody-drug conjugates.

Usage Research use only

Conjugate Biotinylated





Biotinylated Anti-Exatecan antibody(1A6), Rabbit mAb ELISA 0.1µg of BSA-Exatecan protein per well




Figure 1. ELISA plate pre-coated by 1 μ g/mL (100 μ L/well) BSA-Exatecan can Biotinylated Anti-Exatecan antibody(1A6), Rabbit mAb (DME101319B) in a linear range of 16-400 ng/mL.

Email: info@dimabio.com Website: www.dimabio.com

