

PRODUCT INFORMATION

Common Name	JNJ 80202135;JNJ-80202135;JNJ80202135;M281;M-281;M281
Conjugate	Biotinylated
Synonyms	FCRN;FcgammaRn;alpha-chain
Applications	N/A
Endotoxin	Less than 1.0 EU/μg by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	Human IgG1(N297A)-lambda3
Reactivity	Human
Target	FCGRT
Uniprot ID	P55899
Description	Biotinylated Anti-FCGRT(nipocalimab biosimilar) mAb
Delivery	6-8 weeks
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

