

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | BFCR-4350-A, BFCR-4350A, RG-6160, RO-7187797 |
| Conjugate | Biotinylated |
| Synonyms | BXMAS1;CD307;CD307e;FCRH5;IRTA2;PRO820 |
| Applications | ELISA; Flow Cyt |
| Recommended Dilutions | ELISA 1:5000-10000; Flow Cyt 1:100 |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | Humanized |
| IgG type | Human IgG1 - kappa |
| Reactivity | Human |
| Target | FCRL5 |
| Uniprot ID | Q96RD9 |
| Description | Biotinylated Anti-FCRL5(cevostamab biosimilar without anti-CD3 portion) mAb |
| Delivery | 2-3 weeks |
| Storage&Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

