

PRODUCT INFORMATION

| | |
|---|--|
| Common Name | OTSA101, ONC-005 |
| Conjugate | Biotinylated |
| Synonyms | CD350, FZ-10, Fz10, FzE7, hFz10 |
| Applications | ELISA, Flow Cyt |
| Recommended Dilutions | ELISA 1:5000-10000, Flow Cyt 1:100 |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. |
| Host Species | Chimeric |
| IgG type | Human IgG1(E356D,M358L) - kappa |
| Reactivity | Human |
| Target | FZD10 |
| Uniprot ID | Q9ULW2 |
| Description | Biotinylated Anti-FZD10(tabituximab biosimilar) mAb |
| Delivery | 2-3 weeks |
| Storage&Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).Lyophilized antibodies are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

