

PRODUCT INFORMATION

| | |
|------------------------------|---|
| Common Name | PF-06741086;PF06741086;PF-6741086;PF6741086 |
| Conjugate | Biotinylated |
| Synonyms | NEU;NGL;HER2;TKR1;CD340;HER-2;VSCN2;MLN 19;MLN-19;c-ERB2;c-ERB-2;HER-2/neu;p185(erbB2) |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | N/A |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | Her2 |
| Uniprot ID | P04626 |
| Description | Biotinylated Anti-Her2(marstacimab biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage & Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

