

PRODUCT INFORMATION

Common Name AV-299, SCH 900105

Conjugate Biotinylated

Synonyms Hepatopoietin-A, SF, HPTA

Applications ELISA

Recommended

ELISA 1:5000-10000 **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species Humanized

IgG type lgG1 Reactivity Human **Target HGF Uniprot ID** P14210

Biotinylated Anti-HGF(ficlatuzumab biosimilar) **Description**

mAb

Delivery 2-3 weeks

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store

> Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**

