

PRODUCT INFORMATION

Common Name Bimzelx;CDP-4940;CDP4940;UCB-4940;UCB4940

Conjugate Biotinylated

Synonyms IL17;CTLA8;IL-17;ILA17;CTLA-8;IL-17A

Applications N/A
Recommended N/A

Dilutions N/

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

Host Species N/A
IgG type IgG1
Reactivity Human
Target IL17A
Uniprot ID 016552

DescriptionBiotinylated Anti-IL17A(bimekizumab biosimilar)

mAb

Delivery 6-8 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

Storage & Shipping intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

