

PRODUCT INFORMATION

Common Name	S-78989;S78989;S 78989;VPM087;VPM087;VPM-087;VPM87;VPM-87;XMA-005.2;XMA-0052;XMA52;XOMA-052;XOMA052;XOMA52;XOMA-52
Conjugate	Biotinylated
Synonyms	IL-1;IL1F2;IL1beta;IL1-BETA
Applications	N/A
Endotoxin	Less than 1.0 EU/μg by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	IgG2
Reactivity	Human
Target	IL1B
Uniprot ID	P01584
Description	Biotinylated Anti-IL1B(gevokizumab biosimilar) mAb
Delivery	6-8 weeks
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μm) prior to use.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

DIMABIO CONFIDENTIAL

