

PRODUCT INFORMATION

|                              |   |
|------------------------------|---|
| Common Name                  | LEO<br>138559;LEO-138559;LEO138559;ARGX112;ARGX-112;LP0145;LP-0145  |
| Conjugate                    | Biotinylated  |
| Synonyms                     | IL22R;CRF2-9;IL22R1   |
| Applications                 | N/A   |
| Recommended Dilutions        | N/A   |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.  |
| Host Species                 | N/A   |
| IgG type                     | Human IgG1(N297Q,E356D,M358L)-lambda2   |
| Reactivity                   | Human   |
| Target                       | IL22RA1   |
| Uniprot ID                   | Q8N6P7  |
| Description                  | Biotinylated Anti-IL22RA1(temtokibart biosimilar) mAb   |
| Delivery                     | 6-8 weeks   |
| Storage&Shipping             | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background                   | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.   |
| Usage                        | Research use only   |

