

PRODUCT INFORMATION

Common Name DI-17E6, DI17E6, EMD 525797

Conjugate Biotinylated

Synonyms Integrin alpha-V;CD51

Applications Flow Cyt

Recommended

Flow Cyt 1:100 **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species Humanized

IgG type lgG2 Reactivity Human **Target ITGAV Uniprot ID** P06756

Biotinylated Anti-ITGAV(abituzumab biosimilar) **Description**

mAb

Delivery 2-3 weeks

Storage & Shipping

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**

