PRODUCT INFORMATION

Common Name TBC

Conjugate Biotinylated

kidney associated antigen 1;kidney-associated **Synonyms** antigen 1;MGC78738;RU2 antisense gene

protein;RÚ2AS

Applications N/A

Recommended

N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A IgG type IgG1 Reactivity Human **Target** KAAG1 Q9UBP8 **Uniprot ID**

Biotinylated Anti-KAAG1(olintatug biosimilar) mAb **Description**

6-8 weeks **Delivery**

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

temperature.

Research grade biosimilar. Not for use in **Background**

therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

