

## **PRODUCT INFORMATION**

**Common Name** STAR-0215;STAR0215;STAR 0215

Conjugate Biotinylated

PKK;PPK;KLK3;PKKD **Synonyms** 

**Applications** N/A Recommended N/A **Dilutions** 

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

**Host Species** N/A IgG type lgG1 Reactivity Human **Target** KLKB1 **Uniprot ID** P03952

Biotinylated Anti-KLKB1(navenibart biosimilar) **Description** 

mAb

**Delivery** 6-8 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage** 

