

## **PRODUCT INFORMATION**

OMP-59R5 **Common Name** Conjugate Biotinylated

**Synonyms** CADASIL; CADASIL1; CASIL; IMF2; LMNS

**Applications** Flow Cyt

Recommended

Flow Cyt 1:100 **Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

**Host Species** Homo sapiens

IgG type lgG2 Reactivity Human **Target** NOTCH3 **Uniprot ID** Q9UM47

Biotinylated Anti-NOTCH3(tarextumab biosimilar) **Description** 

mAb

**Delivery** 2-3 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage** 

