

PRODUCT INFORMATION

Common Name SCT-I10A Conjugate Biotinylated

Synonyms PD1;PD-1;CD279;SLEB2;hPD-1;hPD-l;hSLE1

Applications N/A Recommended N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A IgG type IgG4 Reactivity Human **Target** PD-1 **Uniprot ID** Q15116

Biotinylated Anti-PD-1(finotonlimab biosimilar) **Description**

mAb

Delivery 6-8 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**

