

## **PRODUCT INFORMATION**

Common Name AK 105;AK105;AK 105

Conjugate Biotinylated

**Synonyms** PD1;PD-1;CD279;SLEB2;hPD-1;hPD-l;hSLE1

Applications N/A
Recommended N/A
Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

Host Species N/A
IgG type IgG1
Reactivity Human
Target PD-1
Uniprot ID 015116

**Description**Biotinylated Anti-PD-1(penpulimab biosimilar)

mAb

**Delivery** 6-8 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

Storage & Shipping intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing)

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

**Usage** Research use only

