

PRODUCT INFORMATION

Common Name ADG104; ADG-104; LDP

Conjugate Biotinylated

Synonyms B7-H;B7H1;PDL1;PD-L1;hPD-L1;PDCD1L1;PDCD1LG1

Applications Recommended

N/A Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for Formulation &

Reconstitution

specific instructions of reconstitution.

Host Species N/A IgG type lgG1 Reactivity Human PD-L1 **Target Uniprot ID** Q9NZQ7

Description Biotinylated Anti-PD-L1(lesabelimab biosimilar) mAb

Delivery 6-8 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use

> > Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping within a month, aliquot and store at -80°C (Avoid

repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans or

animals.

Usage Research use only

