

## PRODUCT INFORMATION

<b>Common Name</b>	PR-1808164;PR1808164;ABBV-0805;ABBV0805;ABBV805;ABBV-805;BAN-0805;BAN0805;BAN805;BAN-805;PD-0805;PD0805;PD-805;PD805
<b>Conjugate</b>	Biotinylated
<b>Synonyms</b>	PD1;NACP;PARK1;PARK4
<b>Applications</b>	N/A
<b>Endotoxin</b>	Less than 1.0 EU/µg by the LAL method. For <1 EU/mg requirements, please contact us for customization.
<b>Recommended Dilutions</b>	N/A
<b>Formulation &amp; Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Host Species</b>	N/A
<b>IgG type</b>	IgG4
<b>Reactivity</b>	Human
<b>Target</b>	SNCA
<b>Uniprot ID</b>	P37840
<b>Description</b>	Biotinylated Anti-SNCA(exidavnemab biosimilar) mAb
<b>Delivery</b>	6-8 weeks
<b>Storage&amp;Shipping</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only

DIMABIO CONFIDENTIAL

