

PRODUCT INFORMATION

Common Name	EB 05;EB05;EB-05;EB 5;EB-5;EB5;EB-07;EB07;EB-07;EB-7;EB7;EB 7;NI-0101;NI 0101;NI0101;NI-101;NI101;NI 101
Conjugate	Biotinylated
Synonyms	TOLL;CD284;TLR-4;ARMD10
Applications	N/A
Endotoxin	Less than 1.0 EU/μg by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	Human IgG1(K214R,N325S,L328F) - Kappa
Reactivity	Human
Target	TLR4
Uniprot ID	O00206
Description	Biotinylated Anti-TLR4(paridiprubart biosimilar) mAb
Delivery	6-8 weeks
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μm) prior to use.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

