

PRODUCT INFORMATION

Clone ID	58A12
Target	TweakR
Synonyms	TNFRSF12A;FGF-inducible 14;FN14;TweakR;CD266
Host Species	Rabbit
Description	Biotinylated Anti-TweakR antibody(58A12), IgG1 Chimeric mAb
Delivery	2-3 weeks
Uniprot ID	Q9NP84
IgG type	Rabbit/Human Fc chimeric IgG1
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1/100
Purification	Purified from cell culture supernatant by affinity chromatography
Endotoxin	Less than 1.0 EU/μg by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μm) prior to use.
Background	TNFRSF12A (TNF Receptor Superfamily Member 12A) is a Protein Coding gene. Diseases associated with TNFRSF12A include Multiple Sclerosis and Stroke;Ischemic. Among its related pathways are ERK Signaling and CDK-mediated phosphorylation and removal of Cdc6
Usage	Research use only
Conjugate	Biotinylated
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scr

